


EasySoaker Hoof Boot Measurement Log

Keep track of your horse's hoof growth.
Questions? Give us a call! 1.800.447.8836.


EasySoaker

Size	Width in inches	Length in inches
Small	3 3/4" - 4 1/2"	4 1/4" - 5"
Medium	4 1/2" - 5 1/4"	5" - 5 3/4"
Large	5 1/4" - 5 3/4"	5 3/4" - 6 1/4"
X-Large	5 3/4" - 6 1/4"	6 1/4" - 6 3/4"
XX-Large	6 1/4" - 6 3/4"	6 3/4" - 7 1/4"
XXX-Large	6 3/4" - 7 1/4"	7 1/4" - 7 3/4"

Size	Width in mm	Length in mm
Small	95-115mm	108-127mm
Medium	115-134mm	127-146mm
Large	134-136mm	146-159mm
X-Large	146-159mm	159-172mm
XX-Large	159-172mm	172-184mm
XXX-Large	172-184mm	184-190mm



Measure the width of the hoof across the bottom at the widest point.



Measure the length of the hoof from the toe to the butters line of the heel. The butters line is the farthest weight bearing point of the heel where the hoof wall ends. Do not include the heel bulbs in the measurement.

HORSE'S NAME:

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			



EasyCare

2300 E. Vistoso Commerce Loop Rd. Tucson, AZ 85755

| 800.447.8836

| 520.297.1900

| www.easycareinc.com